

CHEMBIOCHEM

Supporting Information

© Copyright Wiley-VCH Verlag GmbH & Co. KGaA, 69451 Weinheim, 2006

CHEMBIOCHEM

Supporting Information

for

Tunable DNA Cleavage by Intercalating Peptidoconjugates

Kerry P. Mahon Jr., Marc D. Roy, Jay R. Carreon, Erin G. Prestwich, Jessica L. Rouge,
Stephanie Shin and Shana O. Kelley*

Materials

Solvents were purchased from Fisher and reagents were purchased from Aldrich Chemical Co. (Milwaukee, WI) or Acros Organics (Morris Plains, NJ). Amino acids and Fmoc-Lys(Boc)-Wang resin (75-100 mesh) with a substitution level of 0.70 mmol/g were purchased from Advanced ChemTech (Louisville, KY). Fmoc-protected Rink amide resin (100-200 mesh) was purchased from NovaBiochem (San Diego, CA) with a substitution level of 0.64 mmol/g. Fmoc-Lys(Boc)-OH-Wang resin (75-100 mesh) was purchased from Advanced ChemTech (Louisville, KY) with a substitution level of 0.70 mmol/g. All solvents and reagents were used without further purification. HPLC grade acetonitrile and Millipore water were used for HPLC analysis. Molar extinction coefficient ($\epsilon = 63\,000\text{ M}^{-1}\text{ cm}^{-1}$ in H_2O at 500 nm for TO-conjugates.¹

Analysis

Reversed-phase HPLC (RP-HPLC) was performed with an Agilent 1100 series HPLC using either an analytical column (Varian C18, 5 micron, 250 x 4.6 mm) at a flow rate of 2 mL/min or a preparative column (Hamilton PRP-1, 10 micron, 250 x 10 mm) at

a flow rate of 3.5 mL/min. For analytical HPLC a linear solvent gradient from 5 to 100% B over 60 min was used (solvent A = 0.1% TFA in H₂O; solvent B = 0.1% TFA in CH₃CN). Mass spectral analysis was performed by the Boston College Mass Spectrometry Facility. Samples were analyzed by MALDI-TOF MS on a Micromass ToFSpec 2E using α -cyano-4-hydroxycinnamic acid matrix. UV analysis was performed on an Agilent 8453 Diode Array Spectrophotometer.

General procedure for the synthesis of TO-peptide conjugates

A representative procedure for solid-phase synthesis of TO-conjugate **7** is given below. The term “wash” refers to the following series of rinses: DMF (5 x 1 mL), MeOH (5 x 1 mL) and DCM (5 x 1 mL).

Example synthesis: TO-(d)WK-CONH₂ (7)

Fmoc-Rink amide resin (31.3 mg, 0.0200 mmol, 1.0 equiv) and THF (1 mL) was added to a 1.5 mL Bio-Rad chromatography column. The slurry was agitated for 1 h then the resin was washed.

A 20% piperidine in DMF solution (1 mL) was added to the reaction vessel and agitated for 0.5 h. The resin was then washed.

HBTU (30.3 mg, 0.0800 mmol, 4.0 equiv), Fmoc-D-Lys(Boc)-OH (37.5 mg, 4.0 equiv), dry DMF (1 mL) and *i*-Pr₂NEt (27.9 μ L, 0.160 mmol, 8.0 equiv) was added to the reaction vessel. The reaction was agitated for 3 h. The resin was then washed.

A 20% piperidine in DMF solution (1 mL) was added to the reaction vessel and agitated for 0.5 h. The resin was then washed.

HBTU (30.3 mg, 4.0 equiv), Fmoc-D-Trp(Boc)-OH (42.1 mg, 4.0 equiv), dry DMF (1 mL) and *i*-Pr₂NEt (27.9 μ L, 8.0 equiv) was added to the reaction vessel. The reaction was agitated for 3 h. The resin was then washed.


A 50% piperidine in DMF solution (1 mL) was added to the reaction vessel and agitated for 5 min. The resin was then washed.

HBTU (30.3 mg, 4.0 equiv), TO-COOH (38.8 mg, 4.0 equiv), dry DMF (1 mL) and *i*-Pr₂NEt (27.9 μ L, 8.0 equiv) was added to the reaction vessel, which was then covered with foil. The reaction was agitated for 3 h. The resin was then washed.


10% TFA in DCM (2 x 1 mL) was added to the reaction vessel and the slurry was then transferred to a glass funnel with a fine sinter. The solvent was allowed to drip slowly through the resin bed and was washed with 5% TFA in DCM (2 x 1 mL) and concentrated in vacuo.

Deprotection was achieved by stirring the residue in a 95:5 TFA:TIS solution at room temperature for 0.5 h and concentrating the solution in vacuo. The residue was dissolved in dry MeOH (150 μ L) and transferred to a 15 mL plastic conical tube. Cold Et₂O (10 mL) was added and the resulting crude mixture was centrifuged and the Et₂O was decanted. This procedure was repeated. The red solid was dissolved in 0.1% TFA/H₂O and purified by RP-HPLC.


Characterization of TO-conjugates




TO-GK-COOH (1). MALDI-TOF calculated for [C₃₂H₄₀N₅O₄S⁺] requires *m/z* 590.8 Found 590.8 HPLC retention time 19.6 min




TO-GK-CONH₂ (2). MALDI-TOF calculated for [C₃₂H₄₁N₆O₃S⁺] requires *m/z* 589.8 Found 589.7 HPLC retention time 19.3 min




TO-WK-COOH (**3**). MALDI-TOF calculated for $[C_{41}H_{47}N_6O_4S^+]$ requires m/z 719.9 Found 720.8 HPLC retention time 23.6 min




TO-WK-CONH₂ (**4**). MALDI-TOF calculated for $[C_{41}H_{48}N_7O_3S^+]$ requires m/z 718.9 Found 719.2 HPLC retention time 22.8 min




TO_z-WK-CONH₂ (**5**). MALDI-TOF calculated for $[C_{46}H_{58}N_7O_3S^+]$ requires m/z 789.1 Found 787.9 HPLC retention time 29.1 min




TO_Q-WK-CONH₂ (**6**). MALDI-TOF calculated for [C₄₆H₅₈N₇O₃S⁺] requires *m/z* 789.1
 Found 790.1 HPLC retention time 29.7 min




TO-(d)WK-CONH₂ (**7**). MALDI-TOF calculated for [C₄₁H₄₈N₇O₃S⁺] requires *m/z* 718.9
 Found 719.2 HPLC retention time 22.8 min




TO-WE-CONH₂ (**8**). MALDI-TOF calculated for [C₄₀H₄₃N₆O₅S⁺] requires *m/z* 719.9
 Found 720.2 HPLC retention time 24.5 min




TO-WV-CONH₂ (**9**). MALDI-TOF calculated for [C₄₀H₄₅N₆O₃S⁺] requires *m/z* 689.9
 Found 690.1 HPLC retention time 26.6 min



TO-WQ-CONH₂ (**10**). MALDI-TOF calculated for [C₄₀H₄₄N₇O₄S⁺] requires *m/z* 718.9
 Found 719.1 HPLC retention time 23.6 min



TO-W-CONH₂ (**11**). MALDI-TOF calculated for [C₃₅H₃₆N₅O₂S⁺] requires *m/z* 590.8
 Found 590.8 HPLC retention time 25.7 min



TO-KW-CONH₂ (**12**). MALDI-TOF calculated for [C₄₁H₄₈N₇O₃S⁺] requires *m/z* 718.9
Found 719.2 HPLC retention time 22.9 min

-
- (1) R. P. Haugland, *Handbook of Fluorescent Probes and Research Products*; 9th ed.; Molecular Probes, Eugene, OR, **2002**.